



Our new cardiotoxicity assay, Cardio *quick*Predict, is now available. The assay indicates the cardiotoxicity potential of drug candidates and other lead compounds, based on changes in human induced pluripotent stem cell-derived cardiomyocyte (iPSC-CM) metabolism and viability.

Cardio *quick*Predict provides an indication of cardiotoxicity and identifies functional or structural cardiotoxicants with a single assay. We plan to offer electrophysiological and non-cardiomyocyte viability endpoints in the future. Stay tuned!

[Continued >](#)

iCell® Cardiomyocytes²

Biomarker Response



Evaluate
Cardiotoxicity
Potential across 8-
concentrations

Cardiotoxicity
Prediction

iCell® Cardiomyocytes²

Cell Viability



Non-
Cardiotoxic

Cardiotoxic

Functional

Structural


Future Additions: Non-Cardiomyocyte Cell Viability, Impedance & Electrophysiological Endpoints

DEVELOPMENT PHASE TEST SET PERFORMANCE

Training Set Compounds (N=60)	
POSITIVES (N=40)	NEGATIVES (N=20)
STRUCTURAL Dexfenfluramine Doxorubicin Rofecoxib	Acyclovir Amoxicillin Aspartame Erlotinib
FUNCTIONAL Astemizole Dofetilide Vandetanib	Methapyrilene Ranitidine Sorbitol
STRUCTURAL & FUNCTIONAL Clozapine Isoproterenol Sunitinib	


Test Set Predictions with Biomarker/Viability-Based Model at or within 10-fold of the Therapeutic C _{max}					
True Negatives			True Positives		
Treatment	Dose (µM)	Prediction	Treatment	Dose (µM)	Prediction
Acetylsalicylic Acid	10	Negative	Bepidil	3.3	Positive
Cetirizine	0.8	Negative	Busulfan	150	Positive
Loratadine	0.02	Negative	Levomethadyl Acetate	6	Positive
Sildenafil	1	Negative	Nilotinib	3	Positive
Sucrose	1.8	Negative	Rosiglitazone	1.7	Negative
Xylitol	0.5	Negative	Sotalol	45	Positive

Preliminary model was 92% predictive of cardiotoxicity (100% specificity, 83% sensitivity).



Cell culture

The cardiotoxicity assay is performed using human iPSC-derived cardiomyocytes.



Ultra Performance Liquid Chromatography-Mass Spectrometry

Samples are analyzed using UPLC-HRMS methods that have been optimized for the cardiotoxicity biomarkers of interest, yielding selective and reproducible biomarker measurement.

Data Analysis & Reporting

- SOP-driven analysis
- LIMS-controlled data analysis pipeline
- Identification of critical exposure where cellular metabolism is altered
- Uniform reporting for rapid turnaround
- Custom reporting available

Quality

From start to finish, Stemina has a well-defined, quality program that ensures the integrity of our data.

Experience Counts

Our team has extensive experience in screening a wide variety of proprietary compounds including pharmaceuticals, agri-chemicals, tobacco products, consumer products, and cosmetic ingredients. Stemina was founded in 2006; its state-of-the-art facilities are located in the United States.

Extending Our Global Reach

Stemina has partnered with CiToxLAB, which has facilities in Canada, France, Denmark, and Hungary, to provide worldwide service.

EPA ToxCast™ Contractor